


Set RoSA® IH INSTRUMENTS Removal

Rotation-stable Screw Anchor for extramedullary treatment

non-sterile packaging Set no. 19.610.00



 Internal hexagon
Set no. 19.610.00
non-sterile packaging

Instruments in sets, supplied **non-sterile** Set no. 19.610.00

Item no.	Item no. sterile	Item description	Quant.
Instruments			
2.9406.35/260	—	Hex screwdriver with handle for screws Ø 4.5 mm – 7.0 mm; wrench width 3.5 mm; length 260 mm	1
10.269.80	—	Holding sleeve Ø 12 mm x 130 mm for RoSA® connecting screw	1
10.269.63	—	Extraction rod for RoSA® blade	1
10.269.631	—	Mass for extraction rod for RoSA® blade	1
10.269.632	—	Screw; internal hexagon 3.5 mm; for mass for extraction rod RoSA® blade	1
10.269.41	—	Screwdriver with T-handle for removal RoSA® screw; length 285 mm	1
10.269.50	—	Threaded rod M4 x 305.5 mm for RoSA®; stainless steel	1
10.269.51	—	Screwing aid for rod M4 for RoSA®	1
10.269.72	—	Compression nut; M4 for RoSA®	1
Container			
19.609.00	—	Perforated autoclavable container with inset for implant removal RoSA®	1

Applicable with:

19.606.00	RoSA® IH T SET
19.608.00	RoSA® IH INSTRUMENTS
19.612.00	RoSA® IH T Trochanter stabilization plate SET