

Cannulated screw

titanium

non-sterile packaging / sterile packaging



Art.-Nr. 3.448.30



Art.-Nr. 3.450.30



Art.-Nr. 3.504.30



Art.-Nr. 3.602.30



Art.-Nr. 3.440.30



Art.-Nr. 3.445.30



Art.-Nr. 3.434.30



































Art.-Nr. 3.435.30

| Item no. | EAN | Item description | Item no. sterile | EAN sterile |
|--|-----|---|------------------|-------------|
| Cannulated screw 2.7 cann. 1.5 partial thread st sd | | | | |
| 3.448.10 | | Cannulated screw; Ø 2.7 mm x 10 mm; thread length 5 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.10 | |
| 3.448.12 | | Cannulated screw; Ø 2.7 mm x 12 mm; thread length 5 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.12 | |
| 3.448.14 | | Cannulated screw; Ø 2.7 mm x 14 mm; thread length 5 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.14 | |
| 3.448.16 | | Cannulated screw; Ø 2.7 mm x 16 mm; thread length 6 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.16 | |
| 3.448.18 | | Cannulated screw; Ø 2.7 mm x 18 mm; thread length 7 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.18 | |
| 3.448.20 | | Cannulated screw; Ø 2.7 mm x 20 mm; thread length 8 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.20 | |
| 3.448.22 | | Cannulated screw; Ø 2.7 mm x 22 mm; thread length 9 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.22 | |
| 3.448.24 | | Cannulated screw; Ø 2.7 mm x 24 mm; thread length 10 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.24 | |
| 3.448.26 | | Cannulated screw; Ø 2.7 mm x 26 mm; thread length 10 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.26 | |
| 3.448.28 | | Cannulated screw; Ø 2.7 mm x 28 mm; thread length 12 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.28 | |
| 3.448.30 | | Cannulated screw; Ø 2.7 mm x 30 mm; thread length 12 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.30 | |
| 3.448.32 | | Cannulated screw; Ø 2.7 mm x 32 mm; thread length 12 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.32 | |
| 3.448.34 | | Cannulated screw; Ø 2.7 mm x 34 mm; thread length 12 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.34 | |

Cannulated screw

titanium

































non-sterile packaging / sterile packaging

| Item no. | EAN | Item description | Item no. sterile | EAN sterile |
|----------|---|---|------------------|---|
| 3.448.35 |  | Cannulated screw; Ø 2.7 mm x 35 mm; thread length 15 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.35 |  |
| 3.448.36 |  | Cannulated screw; Ø 2.7 mm x 36 mm; thread length 15 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.36 |  |
| 3.448.38 |  | Cannulated screw; Ø 2.7 mm x 38 mm; thread length 15 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.38 |  |
| 3.448.40 |  | Cannulated screw; Ø 2.7 mm x 40 mm; thread length 15 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.40 |  |
| 3.448.42 |  | Cannulated screw; Ø 2.7 mm x 42 mm; thread length 15 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.42 |  |
| 3.448.44 |  | Cannulated screw; Ø 2.7 mm x 44 mm; thread length 15 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.44 |  |
| 3.448.45 |  | Cannulated screw; Ø 2.7 mm x 45 mm; thread length 15 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.45 |  |
| 3.448.46 |  | Cannulated screw; Ø 2.7 mm x 46 mm; thread length 15 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.46 |  |
| 3.448.48 |  | Cannulated screw; Ø 2.7 mm x 48 mm; thread length 15 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.48 |  |
| 3.448.50 |  | Cannulated screw; Ø 2.7 mm x 50 mm; thread length 15 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.50 |  |
| 3.448.55 |  | Cannulated screw; Ø 2.7 mm x 55 mm; thread length 20 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.55 |  |
| 3.448.60 |  | Cannulated screw; Ø 2.7 mm x 60 mm; thread length 20 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.60 |  |
| 3.448.65 |  | Cannulated screw; Ø 2.7 mm x 65 mm; thread length 20 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.65 |  |
| 3.448.70 |  | Cannulated screw; Ø 2.7 mm x 70 mm; thread length 20 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.70 |  |
| 3.448.75 |  | Cannulated screw; Ø 2.7 mm x 75 mm; thread length 20 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.75 |  |
| 3.448.80 |  | Cannulated screw; Ø 2.7 mm x 80 mm; thread length 20 mm; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.448.80 |  |

Cannulated screw

titanium



































non-sterile packaging / sterile packaging

| Item no. | EAN | Item description | Item no. sterile | EAN sterile |
|--|---|--|------------------|---|
| Cannulated screw 2.7 cann. 1.5 partial thread st sd | | | | |
| 3.450.06 |  | Cannulated screw; Ø 2.7 mm x 6 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.06 |  |
| 3.450.08 |  | Cannulated screw; Ø 2.7 mm x 8 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.08 |  |
| 3.450.10 |  | Cannulated screw; Ø 2.7 mm x 10 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.10 |  |
| 3.450.12 |  | Cannulated screw; Ø 2.7 mm x 12 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.12 |  |
| 3.450.14 |  | Cannulated screw; Ø 2.7 mm x 14 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.14 |  |
| 3.450.16 |  | Cannulated screw; Ø 2.7 mm x 16 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.16 |  |
| 3.450.18 |  | Cannulated screw; Ø 2.7 mm x 18 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.18 |  |
| 3.450.20 |  | Cannulated screw; Ø 2.7 mm x 20 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.20 |  |
| 3.450.22 |  | Cannulated screw; Ø 2.7 mm x 22 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.22 |  |
| 3.450.24 |  | Cannulated screw; Ø 2.7 mm x 24 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.24 |  |
| 3.450.26 |  | Cannulated screw; Ø 2.7 mm x 26 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.26 |  |
| 3.450.28 |  | Cannulated screw; Ø 2.7 mm x 28 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.28 |  |
| 3.450.30 |  | Cannulated screw; Ø 2.7 mm x 30 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.30 |  |
| 3.450.32 |  | Cannulated screw; Ø 2.7 mm x 32 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.32 |  |
| 3.450.34 |  | Cannulated screw; Ø 2.7 mm x 34 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.34 |  |
| 3.450.35 |  | Cannulated screw; Ø 2.7 mm x 35 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.35 |  |

Cannulated screw

titanium

































non-sterile packaging / sterile packaging

| Item no. | EAN | Item description | Item no. sterile | EAN sterile |
|----------|---|---|------------------|---|
| 3.450.36 |  | Cannulated screw; Ø 2.7 mm x 36 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.36 |  |
| 3.450.38 |  | Cannulated screw; Ø 2.7 mm x 38 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.38 |  |
| 3.450.40 |  | Cannulated screw; Ø 2.7 mm x 40 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.40 |  |
| 3.450.42 |  | Cannulated screw; Ø 2.7 mm x 42 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.42 |  |
| 3.450.44 |  | Cannulated screw; Ø 2.7 mm x 44 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.44 |  |
| 3.450.45 |  | Cannulated screw; Ø 2.7 mm x 45 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.45 |  |
| 3.450.46 |  | Cannulated screw; Ø 2.7 mm x 46 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.46 |  |
| 3.450.48 |  | Cannulated screw; Ø 2.7 mm x 48 mm; standard; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.48 |  |
| 3.450.50 |  | Cannulated screw; Ø 2.7 mm x 50 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.50 |  |
| 3.450.55 |  | Cannulated screw; Ø 2.7 mm x 55 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.55 |  |
| 3.450.60 |  | Cannulated screw; Ø 2.7 mm x 60 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.60 |  |
| 3.450.65 |  | Cannulated screw; Ø 2.7 mm x 65 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.65 |  |
| 3.450.70 |  | Cannulated screw; Ø 2.7 mm x 70 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.70 |  |
| 3.450.75 |  | Cannulated screw; Ø 2.7 mm x 75 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.75 |  |
| 3.450.80 |  | Cannulated screw; Ø 2.7 mm x 80 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.80 |  |
| 3.450.65 |  | Cannulated screw; Ø 2.7 mm x 65 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.65 |  |
| 3.450.70 |  | Cannulated screw; Ø 2.7 mm x 70 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB [®] ; titanium | S3.450.70 |  |

Cannulated screw

titanium

































non-sterile packaging / sterile packaging

| Item no. | EAN | Item description | Item no. sterile | EAN sterile |
|---|---|--|------------------|---|
| 3.450.75 |  | Cannulated screw; Ø 2.7 mm x 75 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.75 |  |
| 3.450.80 |  | Cannulated screw; Ø 2.7 mm x 80 mm; fully threaded; internal Ø 1.5 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.450.80 |  |
| Cannulated screw 3.5 cann. 1.3 partial thread st bt sd | | | | |
| 3.504.10 |  | Cannulated screw; Ø 3.5 mm x 10 mm; thread length 5 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.10 |  |
| 3.504.12 |  | Cannulated screw; Ø 3.5 mm x 12 mm; thread length 5 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.12 |  |
| 3.504.14 |  | Cannulated screw; Ø 3.5 mm x 14 mm; thread length 5 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.14 |  |
| 3.504.16 |  | Cannulated screw; Ø 3.5 mm x 16 mm; thread length 6 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.16 |  |
| 3.504.18 |  | Cannulated screw; Ø 3.5 mm x 18 mm; thread length 7 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.18 |  |
| 3.504.20 |  | Cannulated screw; Ø 3.5 mm x 20 mm; thread length 9 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.20 |  |
| 3.504.22 |  | Cannulated screw; Ø 3.5 mm x 22 mm; thread length 10 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.22 |  |
| 3.504.24 |  | Cannulated screw; Ø 3.5 mm x 24 mm; thread length 12 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.24 |  |
| 3.504.26 |  | Cannulated screw; Ø 3.5 mm x 26 mm; thread length 14 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.26 |  |
| 3.504.28 |  | Cannulated screw; Ø 3.5 mm x 28 mm; thread length 14 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.28 |  |
| 3.504.30 |  | Cannulated screw; Ø 3.5 mm x 30 mm; thread length 14 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.30 |  |
| 3.504.32 |  | Cannulated screw; Ø 3.5 mm x 32 mm; thread length 14 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.32 |  |
| 3.504.34 |  | Cannulated screw; Ø 3.5 mm x 34 mm; thread length 14 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.34 |  |
| 3.504.36 |  | Cannulated screw; Ø 3.5 mm x 36 mm; thread length 14 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.36 |  |

Cannulated screw

titanium


non-sterile packaging / sterile packaging

| Item no. | EAN | Item description | Item no. sterile | EAN sterile |
|----------|---|--|------------------|---|
| 3.504.38 |  | Cannulated screw; Ø 3.5 mm x 38 mm; thread length 14 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.38 |  |
| 3.504.40 |  | Cannulated screw; Ø 3.5 mm x 40 mm; thread length 14 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.40 |  |
| 3.504.42 |  | Cannulated screw; Ø 3.5 mm x 42 mm; thread length 15 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.42 |  |
| 3.504.44 |  | Cannulated screw; Ø 3.5 mm x 44 mm; thread length 15 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.44 |  |
| 3.504.46 |  | Cannulated screw; Ø 3.5 mm x 46 mm; thread length 15 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.46 |  |
| 3.504.48 |  | Cannulated screw; Ø 3.5 mm x 48 mm; thread length 15 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.48 |  |
| 3.504.50 |  | Cannulated screw; Ø 3.5 mm x 50 mm; thread length 15 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.50 |  |
| 3.504.52 |  | Cannulated screw; Ø 3.5 mm x 52 mm; thread length 16 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.52 |  |
| 3.504.54 |  | Cannulated screw; Ø 3.5 mm x 54 mm; thread length 16 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.54 |  |
| 3.504.56 |  | Cannulated screw; Ø 3.5 mm x 56 mm; thread length 16 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.56 |  |
| 3.504.58 |  | Cannulated screw; Ø 3.5 mm x 58 mm; thread length 16 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.58 |  |
| 3.504.60 |  | Cannulated screw; Ø 3.5 mm x 58 mm; thread length 16 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.60 |  |
| 3.504.62 |  | Cannulated screw; Ø 3.5 mm x 60 mm; thread length 16 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.62 |  |
| 3.504.64 |  | Cannulated screw; Ø 3.5 mm x 62 mm; thread length 17 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.64 |  |
| 3.504.66 |  | Cannulated screw; Ø 3.5 mm x 66 mm; thread length 17 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.66 |  |
| 3.504.68 |  | Cannulated screw; Ø 3.5 mm x 68 mm; thread length 17 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.68 |  |

Cannulated screw

titanium

































non-sterile packaging / sterile packaging

| Item no. | EAN | Item description | Item no. sterile | EAN sterile |
|---|---|--|------------------|---|
| 3.504.70 |  | Cannulated screw; Ø 3.5 mm x 70 mm; thread length 17 mm; internal Ø 1.3 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.504.70 |  |
| Cannulated screw 3.5 cann. 1.3 st sd | | | | |
| 3.604.10 |  | Cannulated screw; Ø 3.5 mm x 10 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.10 |  |
| 3.604.12 |  | Cannulated screw; Ø 3.5 mm x 12 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.12 |  |
| 3.604.14 |  | Cannulated screw; Ø 3.5 mm x 14 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.14 |  |
| 3.604.16 |  | Cannulated screw; Ø 3.5 mm x 16 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.16 |  |
| 3.604.18 |  | Cannulated screw; Ø 3.5 mm x 18 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.18 |  |
| 3.604.20 |  | Cannulated screw; Ø 3.5 mm x 20 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.20 |  |
| 3.604.22 |  | Cannulated screw; Ø 3.5 mm x 22 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.22 |  |
| 3.604.24 |  | Cannulated screw; Ø 3.5 mm x 24 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.24 |  |
| 3.604.26 |  | Cannulated screw; Ø 3.5 mm x 26 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.26 |  |
| 3.604.28 |  | Cannulated screw; Ø 3.5 mm x 28 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.28 |  |
| 3.604.30 |  | Cannulated screw; Ø 3.5 mm x 30 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.30 |  |
| 3.604.32 |  | Cannulated screw; Ø 3.5 mm x 32 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.32 |  |
| 3.604.34 |  | Cannulated screw; Ø 3.5 mm x 34 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.34 |  |
| 3.604.36 |  | Cannulated screw; Ø 3.5 mm x 36 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.36 |  |
| 3.604.38 |  | Cannulated screw; Ø 3.5 mm x 38 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.38 |  |

Cannulated screw

titanium

non-sterile packaging / sterile packaging

| Item no. | EAN | Item description | Item no. sterile | EAN sterile |
|----------|---|---|------------------|---|
| 3.604.40 |  | Cannulated screw; Ø 3.5 mm x 40 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.40 |  |
| 3.604.42 |  | Cannulated screw; Ø 3.5 mm x 42 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.42 |  |
| 3.604.44 |  | Cannulated screw; Ø 3.5 mm x 44 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.44 |  |
| 3.604.46 |  | Cannulated screw; Ø 3.5 mm x 46 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.46 |  |
| 3.604.48 |  | Cannulated screw; Ø 3.5 mm x 48 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.48 |  |
| 3.604.50 |  | Cannulated screw; Ø 3.5 mm x 50 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.50 |  |
| 3.604.52 |  | Cannulated screw; Ø 3.5 mm x 52 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.52 |  |
| 3.604.54 |  | Cannulated screw; Ø 3.5 mm x 54 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.54 |  |
| 3.604.56 |  | Cannulated screw; Ø 3.5 mm x 56 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.56 |  |
| 3.604.58 |  | Cannulated screw; Ø 3.5 mm x 58 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.58 |  |
| 3.604.60 |  | Cannulated screw; Ø 3.5 mm x 60 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.60 |  |
| 3.604.62 |  | Cannulated screw; Ø 3.5 mm x 62 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.62 |  |
| 3.604.64 |  | Cannulated screw; Ø 3.5 mm x 64 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.64 |  |
| 3.604.66 |  | Cannulated screw; Ø 3.5 mm x 66 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.66 |  |
| 3.604.68 |  | Cannulated screw; Ø 3.5 mm x 68 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.68 |  |
| 3.604.70 |  | Cannulated screw; Ø 3.5 mm x 70 mm; fully threaded; internal Ø 1.3 mm; self drilling; self tapping; internal hexagon 2.5 mm; titanium | S3.604.70 |  |

Cannulated screw 4.0 cann. 1.8 partial thread st bt sd



Königsee Implantate GmbH
OT Aschau | Am Sand 4
07426 Allendorf | Germany

Fon +49 36738 498-560
Fax +49 36738 498-569

Mail export@koenigsee-implantate.de



































For further information please visit:
www.koenigsee-implantate.de

AnWa | Edition 05 | 2025-10
All pictured products are CE compliant.

Cannulated screw

titanium

non-sterile packaging / sterile packaging

| Item no. | EAN | Item description | Item no. sterile | EAN sterile |
|----------|---|---|------------------|---|
| 3.440.10 |  | Cannulated screw; Ø 4.0 mm x 10 mm; thread length 5 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.10 |  |
| 3.440.12 |  | Cannulated screw; Ø 4.0 mm x 12 mm; thread length 5 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.12 |  |
| 3.440.14 |  | Cannulated screw; Ø 4.0 mm x 14 mm; thread length 5 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.14 |  |
| 3.440.16 |  | Cannulated screw; Ø 4.0 mm x 16 mm; thread length 6 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.16 |  |
| 3.440.18 |  | Cannulated screw; Ø 4.0 mm x 18 mm; thread length 7 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.18 |  |
| 3.440.20 |  | Cannulated screw; Ø 4.0 mm x 20 mm; thread length 8 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.20 |  |
| 3.440.22 |  | Cannulated screw; Ø 4.0 mm x 22 mm; thread length 9 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.22 |  |
| 3.440.24 |  | Cannulated screw; Ø 4.0 mm x 24 mm; thread length 10 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.24 |  |
| 3.440.26 |  | Cannulated screw; Ø 4.0 mm x 26 mm; thread length 12 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.26 |  |
| 3.440.28 |  | Cannulated screw; Ø 4.0 mm x 28 mm; thread length 14 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.28 |  |
| 3.440.30 |  | Cannulated screw; Ø 4.0 mm x 30 mm; thread length 14 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.30 |  |
| 3.440.32 |  | Cannulated screw; Ø 4.0 mm x 32 mm; thread length 14 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.32 |  |
| 3.440.34 |  | Cannulated screw; Ø 4.0 mm x 34 mm; thread length 14 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.34 |  |
| 3.440.35 |  | Cannulated screw; Ø 4.0 mm x 35 mm; thread length 14 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.35 |  |
| 3.440.36 |  | Cannulated screw; Ø 4.0 mm x 36 mm; thread length 14 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.36 |  |
| 3.440.38 |  | Cannulated screw; Ø 4.0 mm x 38 mm; thread length 14 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.38 |  |
| 3.440.40 |  | Cannulated screw; Ø 4.0 mm x 40 mm; thread length 14 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.40 |  |









Cannulated screw

titanium

non-sterile packaging / sterile packaging

| Item no. | EAN | Item description | Item no. sterile | EAN sterile |
|----------|---|---|------------------|---|
| 3.440.42 |  | Cannulated screw; Ø 4.0 mm x 42 mm; thread length 15 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.42 |  |
| 3.440.44 |  | Cannulated screw; Ø 4.0 mm x 44 mm; thread length 15 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.44 |  |
| 3.440.45 |  | Cannulated screw; Ø 4.0 mm x 45 mm; thread length 15 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.45 |  |
| 3.440.46 |  | Cannulated screw; Ø 4.0 mm x 46 mm; thread length 15 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.46 |  |
| 3.440.48 |  | Cannulated screw; Ø 4.0 mm x 48 mm; thread length 15 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.48 |  |
| 3.440.50 |  | Cannulated screw; Ø 4.0 mm x 50 mm; thread length 15 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.50 |  |
| 3.440.55 |  | Cannulated screw; Ø 4.0 mm x 55 mm; thread length 16 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.55 |  |
| 3.440.60 |  | Cannulated screw; Ø 4.0 mm x 60 mm; thread length 16 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.60 |  |
| 3.440.65 |  | Cannulated screw; Ø 4.0 mm x 65 mm; thread length 17 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.65 |  |
| 3.440.70 |  | Cannulated screw; Ø 4.0 mm x 70 mm; thread length 17 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.70 |  |
| 3.440.75 |  | Cannulated screw; Ø 4.0 mm x 75 mm; thread length 17 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.75 |  |
| 3.440.80 |  | Cannulated screw; Ø 4.0 mm x 80 mm; thread length 18 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.440.80 |  |



































Cannulated screw 4.0 cann. 1.8 st sd

| | | | | |
|----------|---|--|-----------|---|
| 3.445.16 |  | Cannulated screw; Ø 4.0 mm x 16 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.16 |  |
| 3.445.18 |  | Cannulated screw; Ø 4.0 mm x 18 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.18 |  |
| 3.445.20 |  | Cannulated screw; Ø 4.0 mm x 20 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.18 |  |
| 3.445.22 |  | Cannulated screw; Ø 4.0 mm x 22 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.22 |  |

Cannulated screw

titanium









non-sterile packaging / sterile packaging

| Art.-Nr. | EAN | Beschreibung | Art.-Nr. steril | EAN steril |
|----------|---|--|-----------------|---|
| 3.445.24 |  | Cannulated screw; Ø 4.0 mm x 24 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.24 |  |
| 3.445.26 |  | Cannulated screw; Ø 4.0 mm x 26 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.26 |  |
| 3.445.28 |  | Cannulated screw; Ø 4.0 mm x 28 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.28 |  |
| 3.445.30 |  | Cannulated screw; Ø 4.0 mm x 30 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.30 |  |
| 3.445.32 |  | Cannulated screw; Ø 4.0 mm x 32 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.32 |  |
| 3.445.34 |  | Cannulated screw; Ø 4.0 mm x 34 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.34 |  |
| 3.445.35 |  | Cannulated screw; Ø 4.0 mm x 35 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.35 |  |
| 3.445.36 |  | Cannulated screw; Ø 4.0 mm x 36 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.36 |  |
| 3.445.38 |  | Cannulated screw; Ø 4.0 mm x 38 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.38 |  |
| 3.445.40 |  | Cannulated screw; Ø 4.0 mm x 40 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.40 |  |
| 3.445.42 |  | Cannulated screw; Ø 4.0 mm x 42 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.42 |  |
| 3.445.44 |  | Cannulated screw; Ø 4.0 mm x 44 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.44 |  |
| 3.445.45 |  | Cannulated screw; Ø 4.0 mm x 45 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.45 |  |
| 3.445.46 |  | Cannulated screw; Ø 4.0 mm x 46 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.46 |  |
| 3.445.48 |  | Cannulated screw; Ø 4.0 mm x 48 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.48 |  |
| 3.445.50 |  | Cannulated screw; Ø 4.0 mm x 50 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.50 |  |
| 3.445.55 |  | Cannulated screw; Ø 4.0 mm x 55 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.55 |  |

Cannulated screw

titanium































non-sterile packaging / sterile packaging

| Art.-Nr. | EAN | Beschreibung | Art.-Nr. steril | EAN steril |
|---|---|--|-----------------|---|
| 3.445.60 |  | Cannulated screw; Ø 4.0 mm x 60 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.60 |  |
| 3.445.65 |  | Cannulated screw; Ø 4.0 mm x 65 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.65 |  |
| 3.445.70 |  | Cannulated screw; Ø 4.0 mm x 70 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 2.5 mm; TiOB®; titanium | S3.445.70 |  |
| Cannulated screw 4.5 cann. 1.8 partial thread st bt sd | | | | |
| 3.434.20 |  | Cannulated screw; Ø 4.5 mm x 20 mm; thread length 7 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.20 |  |
| 3.434.22 |  | Cannulated screw; Ø 4.5 mm x 22 mm; thread length 7 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.22 |  |
| 3.434.24 |  | Cannulated screw; Ø 4.5 mm x 24 mm; thread length 8 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.24 |  |
| 3.434.26 |  | Cannulated screw; Ø 4.5 mm x 26 mm; thread length 9 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.26 |  |
| 3.434.28 |  | Cannulated screw; Ø 4.5 mm x 28 mm; thread length 9 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.28 |  |
| 3.434.30 |  | Cannulated screw; Ø 4.5 mm x 30 mm; thread length 9 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.30 |  |
| 3.434.32 |  | Cannulated screw; Ø 4.5 mm x 32 mm; thread length 9 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.32 |  |
| 3.434.34 |  | Cannulated screw; Ø 4.5 mm x 34 mm; thread length 9 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.34 |  |
| 3.434.36 |  | Cannulated screw; Ø 4.5 mm x 36 mm; thread length 9 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.36 |  |
| 3.434.38 |  | Cannulated screw; Ø 4.5 mm x 38 mm; thread length 9 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.38 |  |
| 3.434.40 |  | Cannulated screw; Ø 4.5 mm x 40 mm; thread length 9 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.40 |  |
| 3.434.42 |  | Cannulated screw; Ø 4.5 mm x 42 mm; thread length 14 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.42 |  |
| 3.434.44 |  | Cannulated screw; Ø 4.5 mm x 44 mm; thread length 15 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.44 |  |

Cannulated screw

titanium



































non-sterile packaging / sterile packaging

| Art.-Nr. | EAN | Beschreibung | Art.-Nr. steril | EAN steril |
|---|---|--|-----------------|---|
| 3.434.46 |  | Cannulated screw; Ø 4.5 mm x 46 mm; thread length 15 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.46 |  |
| 3.434.48 |  | Cannulated screw; Ø 4.5 mm x 48 mm; thread length 16 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.48 |  |
| Cannulated screw 4.5 cann. 1.8 partial thread st bt sd | | | | |
| 3.434.50 |  | Cannulated screw; Ø 4.5 mm x 50 mm; thread length 16 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.50 |  |
| 3.434.52 |  | Cannulated screw; Ø 4.5 mm x 52 mm; thread length 17 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.52 |  |
| 3.434.54 |  | Cannulated screw; Ø 4.5 mm x 54 mm; thread length 19 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.54 |  |
| 3.434.56 |  | Cannulated screw; Ø 4.5 mm x 56 mm; thread length 19 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.56 |  |
| 3.434.58 |  | Cannulated screw; Ø 4.5 mm x 58 mm; thread length 20 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.58 |  |
| 3.434.60 |  | Cannulated screw; Ø 4.5 mm x 60 mm; thread length 20 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.60 |  |
| 3.434.62 |  | Cannulated screw; Ø 4.5 mm x 62 mm; thread length 21 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.62 |  |
| 3.434.64 |  | Cannulated screw; Ø 4.5 mm x 64 mm; thread length 21 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.64 |  |
| 3.434.66 |  | Cannulated screw; Ø 4.5 mm x 66 mm; thread length 22 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.66 |  |
| 3.434.68 |  | Cannulated screw; Ø 4.5 mm x 68 mm; thread length 23 mm; internal Ø 1.8 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.434.68 |  |
| 3.434.72 |  | Cannulated screw; Ø 4.5 mm x 72 mm; thread length 24 mm; internal Ø 1.8 mm; self drilling; back tapping; self tapping; internal hexagon 3.5 mm; titanium | S3.434.72 |  |
| Cannulated screw 4.5 cann. 1.8 st sd | | | | |
| 3.439.20 |  | Cannulated screw; Ø 4.5 mm x 20 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.20 |  |
| 3.439.22 |  | Cannulated screw; Ø 4.5 mm x 22 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.22 |  |

Cannulated screw

titanium







non-sterile packaging / sterile packaging

| Art.-Nr. | EAN | Beschreibung | Art.-Nr. steril | EAN steril |
|----------|---|---|-----------------|---|
| 3.439.24 |  | Cannulated screw; Ø 4.5 mm x 24 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.24 |  |
| 3.439.26 |  | Cannulated screw; Ø 4.5 mm x 26 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.26 |  |
| 3.439.28 |  | Cannulated screw; Ø 4.5 mm x 28 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.28 |  |
| 3.439.30 |  | Cannulated screw; Ø 4.5 mm x 30 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.30 |  |
| 3.439.32 |  | Cannulated screw; Ø 4.5 mm x 32 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.32 |  |
| 3.439.34 |  | Cannulated screw; Ø 4.5 mm x 34 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.34 |  |
| 3.439.36 |  | Cannulated screw; Ø 4.5 mm x 36 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.36 |  |
| 3.439.38 |  | Cannulated screw; Ø 4.5 mm x 36 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.38 |  |
| 3.439.40 |  | Cannulated screw; Ø 4.5 mm x 38 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.40 |  |
| 3.439.42 |  | Cannulated screw; Ø 4.5 mm x 40 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.42 |  |
| 3.439.44 |  | Cannulated screw; Ø 4.5 mm x 42 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.44 |  |
| 3.439.46 |  | Cannulated screw; Ø 4.5 mm x 44 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.46 |  |
| 3.439.48 |  | Cannulated screw; Ø 4.5 mm x 46 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.48 |  |
| 3.439.50 |  | Cannulated screw; Ø 4.5 mm x 48 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.50 |  |
| 3.439.52 |  | Cannulated screw; Ø 4.5 mm x 52 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.52 |  |
| 3.439.54 |  | Cannulated screw; Ø 4.5 mm x 54 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.54 |  |
| 3.439.56 |  | Cannulated screw; Ø 4.5 mm x 56 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.56 |  |



















Cannulated screw

titanium

non-sterile packaging / sterile packaging

| Art.-Nr. | EAN | Beschreibung | Art.-Nr. steril | EAN steril |
|----------|--|---|-----------------|---|
| 3.439.58 |  | Cannulated screw; Ø 4.5 mm x 58 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.58 |  |
| 3.439.60 |  | Cannulated screw; Ø 4.5 mm x 60 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.60 |  |
| 3.439.62 |  | Cannulated screw; Ø 4.5 mm x 62 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.62 |  |
| 3.439.64 |  | Cannulated screw; Ø 4.5 mm x 64 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.64 |  |
| 3.439.66 |  | Cannulated screw; Ø 4.5 mm x 66 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.66 |  |
| 3.439.68 |  | Cannulated screw; Ø 4.5 mm x 68 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.68 |  |
| 3.439.72 |  | Cannulated screw; Ø 4.5 mm x 72 mm; fully threaded; internal Ø 1.8 mm; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.439.72 |  |











Cannulated screw 6.5 cann. 2.2 thread 16 st bt sd

| | | | | |
|----------|---|---|-----------|---|
| 3.520.30 |  | Cannulated screw; Ø 6.5 mm x 30 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.30 |  |
| 3.520.35 |  | Cannulated screw; Ø 6.5 mm x 35 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.35 |  |
| 3.520.40 |  | Cannulated screw; Ø 6.5 mm x 40 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.40 |  |
| 3.520.45 |  | Cannulated screw; Ø 6.5 mm x 45 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.45 |  |
| 3.520.50 |  | Cannulated screw; Ø 6.5 mm x 50 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.50 |  |
| 3.520.55 |  | Cannulated screw; Ø 6.5 mm x 55 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.55 |  |
| 3.520.60 |  | Cannulated screw; Ø 6.5 mm x 60 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.60 |  |
| 3.520.65 |  | Cannulated screw; Ø 6.5 mm x 65 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.65 |  |
| 3.520.70 |  | Cannulated screw; Ø 6.5 mm x 70 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.70 |  |





Cannulated screw

titanium

non-sterile packaging / sterile packaging

| Art.-Nr. | EAN | Beschreibung | Art.-Nr. steril | EAN steril |
|-----------|---|--|-----------------|---|
| 3.520.75 |  | Cannulated screw; Ø 6.5 mm x 75 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.75 |  |
| 3.520.80 |  | Cannulated screw; Ø 6.5 mm x 80 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.80 |  |
| 3.520.85 |  | Cannulated screw; Ø 6.5 mm x 85 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.85 |  |
| 3.520.90 |  | Cannulated screw; Ø 6.5 mm x 90 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.90 |  |
| 3.520.95 |  | Cannulated screw; Ø 6.5 mm x 95 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.95 |  |
| 3.520.100 |  | Cannulated screw; Ø 6.5 mm x 100 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.100 |  |
| 3.520.105 |  | Cannulated screw; Ø 6.5 mm x 105 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.105 |  |
| 3.520.110 |  | Cannulated screw; Ø 6.5 mm x 110 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.110 |  |
| 3.520.115 |  | Cannulated screw; Ø 6.5 mm x 115 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.115 |  |
| 3.520.120 |  | Cannulated screw; Ø 6.5 mm x 120 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.120 |  |
| 3.520.125 |  | Cannulated screw; Ø 6.5 mm x 125 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.125 |  |
| 3.520.130 |  | Cannulated screw; Ø 6.5 mm x 130 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.130 |  |
| 3.520.135 |  | Cannulated screw; Ø 6.5 mm x 135 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.135 |  |
| 3.520.140 |  | Cannulated screw; Ø 6.5 mm x 140 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.520.140 |  |



































Cannulated screw 6.5 cann. 2.2 thread 32 st bt sd

| | | | | |
|----------|---|---|-----------|---|
| 3.523.45 |  | Cannulated screw; Ø 6.5 mm x 45 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.45 |  |
| 3.523.50 |  | Cannulated screw; Ø 6.5 mm x 50 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.50 |  |

Cannulated screw

titanium

non-sterile packaging / sterile packaging

| Art.-Nr. | EAN | Beschreibung | Art.-Nr. steril | EAN steril |
|-----------|---|---|-----------------|---|
| 3.523.55 |  | Cannulated screw; Ø 6.5 mm x 55 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.55 |  |
| 3.523.60 |  | Cannulated screw; Ø 6.5 mm x 60 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.60 |  |
| 3.523.65 |  | Cannulated screw; Ø 6.5 mm x 65 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.65 |  |
| 3.523.70 |  | Cannulated screw; Ø 6.5 mm x 70 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.70 |  |
| 3.523.75 |  | Cannulated screw; Ø 6.5 mm x 75 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.75 |  |
| 3.523.80 |  | Cannulated screw; Ø 6.5 mm x 80 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.80 |  |
| 3.523.85 |  | Cannulated screw; Ø 6.5 mm x 85 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.85 |  |
| 3.523.90 |  | Cannulated screw; Ø 6.5 mm x 90 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.90 |  |
| 3.523.95 |  | Cannulated screw; Ø 6.5 mm x 95 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.95 |  |
| 3.523.100 |  | Cannulated screw; Ø 6.5 mm x 100 mm; standard; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.100 |  |
| 3.523.105 |  | Cannulated screw; Ø 6.5 mm x 105 mm; standard; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.105 |  |
| 3.523.110 |  | Cannulated screw; Ø 6.5 mm x 110 mm; standard; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.110 |  |
| 3.523.115 |  | Cannulated screw; Ø 6.5 mm x 115 mm; standard; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.115 |  |
| 3.523.120 |  | Cannulated screw; Ø 6.5 mm x 120 mm; standard; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.120 |  |
| 3.523.125 |  | Cannulated screw; Ø 6.5 mm x 125 mm; standard; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.125 |  |
| 3.523.130 |  | Cannulated screw; Ø 6.5 mm x 130 mm; standard; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.130 |  |
| 3.523.135 |  | Cannulated screw; Ø 6.5 mm x 135 mm; standard; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.135 |  |

Cannulated screw

titanium






























non-sterile packaging / sterile packaging

| Art.-Nr. | EAN | Beschreibung | Art.-Nr. steril | EAN steril |
|--|---|---|-----------------|---|
| 3.523.140 |  | Cannulated screw; Ø 6.5 mm x 140 mm; standard; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.523.140 |  |
| Cannulated screw 6.5 cann. 2.2 st bt sd | | | | |
| 3.525.45 |  | Cannulated screw; Ø 6.5 mm x 45 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.45 |  |
| 3.525.50 |  | Cannulated screw; Ø 6.5 mm x 50 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.50 |  |
| 3.525.55 |  | Cannulated screw; Ø 6.5 mm x 55 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.55 |  |
| 3.525.60 |  | Cannulated screw; Ø 6.5 mm x 60 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.60 |  |
| 3.525.65 |  | Cannulated screw; Ø 6.5 mm x 65 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.65 |  |
| 3.525.70 |  | Cannulated screw; Ø 6.5 mm x 70 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.70 |  |
| 3.525.75 |  | Cannulated screw; Ø 6.5 mm x 75 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.75 |  |
| 3.525.80 |  | Cannulated screw; Ø 6.5 mm x 80 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.80 |  |
| 3.525.85 |  | Cannulated screw; Ø 6.5 mm x 85 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.85 |  |
| 3.525.90 |  | Cannulated screw; Ø 6.5 mm x 90 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.90 |  |
| 3.525.95 |  | Cannulated screw; Ø 6.5 mm x 95 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.95 |  |
| 3.525.100 |  | Cannulated screw; Ø 6.5 mm x 100 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.100 |  |
| 3.525.105 |  | Cannulated screw; Ø 6.5 mm x 105 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.105 |  |
| 3.525.110 |  | Cannulated screw; Ø 6.5 mm x 110 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.110 |  |
| 3.525.115 |  | Cannulated screw; Ø 6.5 mm x 115 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.115 |  |

Cannulated screw

titanium































non-sterile packaging / sterile packaging

| Art.-Nr. | EAN | Beschreibung | Art.-Nr. steril | EAN steril |
|--|---|---|-----------------|---|
| 3.525.120 |  | Cannulated screw; Ø 6.5 mm x 120 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.120 |  |
| 3.525.125 |  | Cannulated screw; Ø 6.5 mm x 125 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.125 |  |
| 3.525.130 |  | Cannulated screw; Ø 6.5 mm x 130 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.130 |  |
| 3.525.135 |  | Cannulated screw; Ø 6.5 mm x 135 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.135 |  |
| 3.525.140 |  | Cannulated screw; Ø 6.5 mm x 140 mm; fully threaded; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; TiOB®; titanium | S3.525.140 |  |
| Cannulated screw 7.0 cann. 2.2 thread 16 st bt sd | | | | |
| 3.414.30 |  | Cannulated screw; Ø 7.0 mm x 30 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.414.30 |  |
| 3.414.35 |  | Cannulated screw; Ø 7.0 mm x 35 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.414.35 |  |
| 3.414.40 |  | Cannulated screw; Ø 7.0 mm x 40 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.414.40 |  |
| 3.414.45 |  | Cannulated screw; Ø 7.0 mm x 45 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.414.45 |  |
| 3.414.50 |  | Cannulated screw; Ø 7.0 mm x 50 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.414.50 |  |
| 3.414.55 |  | Cannulated screw; Ø 7.0 mm x 55 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.414.55 |  |
| 3.414.60 |  | Cannulated screw; Ø 7.0 mm x 60 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.414.60 |  |
| 3.414.65 |  | Cannulated screw; Ø 7.0 mm x 65 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.414.65 |  |
| 3.414.70 |  | Cannulated screw; Ø 7.0 mm x 70 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.414.70 |  |
| 3.414.75 |  | Cannulated screw; Ø 7.0 mm x 75 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.414.75 |  |
| 3.414.80 |  | Cannulated screw; Ø 7.0 mm x 80 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.414.80 |  |

Cannulated screw

titanium





























non-sterile packaging / sterile packaging

| Art.-Nr. | EAN | Beschreibung | Art.-Nr. steril | EAN steril |
|--|---|--|-----------------|---|
| 3.414.85 |  | Cannulated screw; Ø 7.0 mm x 85 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.414.85 |  |
| 3.414.90 |  | Cannulated screw; Ø 7.0 mm x 90 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.414.90 |  |
| Cannulated screw 7.0 cann. 2.2 thread 16 st bt sd | | | | |
| 3.414.95 |  | Cannulated screw; Ø 7.0 mm x 95 mm; thread length 16 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.414.95 |  |
| 3.414.100 |  | Cannulated screw; Ø 7.0 mm x 100 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self tapping; internal hexagon 3.5 mm; titanium | S3.414.100 |  |
| 3.414.105 |  | Cannulated screw; Ø 7.0 mm x 105 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self tapping; internal hexagon 3.5 mm; titanium | S3.414.105 |  |
| 3.414.110 |  | Cannulated screw; Ø 7.0 mm x 110 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self tapping; internal hexagon 3.5 mm; titanium | S3.414.110 |  |
| 3.414.115 |  | Cannulated screw; Ø 7.0 mm x 115 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self tapping; internal hexagon 3.5 mm; titanium | S3.414.115 |  |
| 3.414.120 |  | Cannulated screw; Ø 7.0 mm x 120 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self tapping; internal hexagon 3.5 mm; titanium | S3.414.120 |  |
| 3.414.125 |  | Cannulated screw; Ø 7.0 mm x 125 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self tapping; internal hexagon 3.5 mm; titanium | S3.414.125 |  |
| 3.414.130 |  | Cannulated screw; Ø 7.0 mm x 130 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self tapping; internal hexagon 3.5 mm; titanium | S3.414.130 |  |
| Cannulated screw 7.0 cann. 2.2 thread 32 st bt sd | | | | |
| 3.424.40 |  | Cannulated screw; Ø 7.0 mm x 40 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.40 |  |
| 3.424.45 |  | Cannulated screw; Ø 7.0 mm x 45 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.45 |  |
| 3.424.50 |  | Cannulated screw; Ø 7.0 mm x 50 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.50 |  |
| 3.424.55 |  | Cannulated screw; Ø 7.0 mm x 55 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.55 |  |
| 3.424.60 |  | Cannulated screw; Ø 7.0 mm x 60 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.60 |  |

Cannulated screw

titanium

non-sterile packaging / sterile packaging

| Art.-Nr. | EAN | Beschreibung | Art.-Nr. steril | EAN steril |
|-----------|---|---|-----------------|---|
| 3.424.65 |  | Cannulated screw; Ø 7.0 mm x 65 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.65 |  |
| 3.424.70 |  | Cannulated screw; Ø 7.0 mm x 70 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.70 |  |
| 3.424.75 |  | Cannulated screw; Ø 7.0 mm x 75 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.75 |  |
| 3.424.80 |  | Cannulated screw; Ø 7.0 mm x 80 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.80 |  |
| 3.424.85 |  | Cannulated screw; Ø 7.0 mm x 85 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.85 |  |
| 3.424.90 |  | Cannulated screw; Ø 7.0 mm x 90 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.90 |  |
| 3.424.95 |  | Cannulated screw; Ø 7.0 mm x 95 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.95 |  |
| 3.424.100 |  | Cannulated screw; Ø 7.0 mm x 100 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.100 |  |
| 3.424.105 |  | Cannulated screw; Ø 7.0 mm x 105 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.105 |  |
| 3.424.110 |  | Cannulated screw; Ø 7.0 mm x 110 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.110 |  |
| 3.424.115 |  | Cannulated screw; Ø 7.0 mm x 115 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.115 |  |
| 3.424.120 |  | Cannulated screw; Ø 7.0 mm x 120 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.120 |  |
| 3.424.125 |  | Cannulated screw; Ø 7.0 mm x 125 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.125 |  |
| 3.424.130 |  | Cannulated screw; Ø 7.0 mm x 130 mm; thread length 32 mm; internal Ø 2.2 mm; back tapping; self drilling; self tapping; internal hexagon 3.5 mm; titanium | S3.424.130 |  |