



# Tibia head plates titanium INSTRUMENTS

for screws  $\varnothing$  4.0 mm

non-sterile packaging Set no. 19.493.40I



 Six lobe  
Set no. 19.493.40I  
non-sterile packaging

Implants in sets, supplied non-sterile Set no. 19.493.40I

Item no.	Item no. sterile	Item description	Quant.
<b>Instruments</b>			
2.904.26	S2.904.26	Drill bit for quick coupling $\varnothing$ 3.0 mm; spiral 60 mm; 2-flute; length 195 mm	1
2.9385.15	—	Screwdriver with handle for six lobe 15; length 215 mm	1
2.953.80	—	Gauge with clasp for screws $\varnothing$ 3.5 mm and $\varnothing$ 4.0 mm; measuring range 90 mm	1
2.977.16	—	Drill guide for angle-stable screwing 4.0 mm; for drills $\varnothing$ 3.0 mm; length 68 mm	2
2.9774.30	—	Sinkhole drill guide for drills $\varnothing$ 3.0 mm/ $\varnothing$ 2.5 mm	1
2.9419.15	—	Screwdriver insert for quick coupling; six lobe 15; length 100 mm	1
<b>Wires</b>			
6.031.18	S6.031.18	Kirschner wire with trocar point and round end; $\varnothing$ 1.8 mm x 150 mm; stainless steel	5
<b>Container</b>			
19.492.40	—	Perforated container with insert for instruments for tibia head plates	1