



Set SF 3.5 tools SL INSTRUMENTS

for screws Ø 4.0/3.5 mm

non-sterile packaging Set no. 19.035.30I



★ Six lobe
Set no. 19.035.30I
Level 1
non-sterile packaging



★ Six lobe
Set no. 19.035.30I
Level 2
non-sterile packaging

Implants in sets, supplied **non-sterile** Set no. 19.035.30I

Item no.	Item no. sterile	Item description	Quant.
Instruments – Level 1			
2.9385.15	–	Screwdriver with handle for six lobe 15; length 215 mm	1
2.904.22	S2.904.22	Drill bit for quick coupling Ø 2.5 mm; spiral 60 mm; 2-flute; length 155 mm; green	1
2.904.25	S2.904.25	Drill bit for quick coupling Ø 3.0 mm; spiral 60 mm; 2-flute; length 155 mm	1
2.977.22	–	Drill guide for angle-stable screwing 4.0; for drills Ø 3.0 mm; double threaded; length 68 mm	2
2.977.21	–	Drill guide with handle variable angle-stable screwing 4.0; for drills Ø 3.0 mm; length 68 mm	1
2.9776.30	–	Sinkhole drill guide with handle; for drills Ø 3.0 mm/2.5 mm	1
Wires			
6.031.16	S6.031.16	Kirschner wire with trocar point and round end; Ø 1.6 mm x 150 mm; stainless steel	5
Instruments – Level 2			
2.977.23	–	Drill guide angle-stable; for drills Ø 2.5 mm; double threaded; length 68 mm; stainless steel	2
10.173.01	–	Sleeve for guide block	2
2.977.09	–	Drill guide with handle for variable angle-stable screwing; for drills Ø 2.5 mm; length 68 mm; small fragment	1
2.953.70	–	Gauge with clasp for screws Ø 3.5 mm and Ø 4.0 mm; measuring range 70 mm	1
10.179.11	–	Aiming block for proximal humerus plate EIGHT ²	1
10.179.16	–	Aiming block for VaWiKo [®] proximal humerus plate	1
10.173.05	–	Screw driver for guide block	1

Set SF 3.5 tools SL INSTRUMENTS


for screws Ø 4.0/3.5 mm

non-sterile packaging Set no. 19.035.30I


Item no.	Item no. sterile	Item description	Quant.
Container			
19.530.10	—	Perforated autoclavable container with inset for instruments-modul small fragment	1

Applicable with:

19.3673.10	Prox. Humerus pl. EIGHT ² vas SL T SET
19.3673.10PS+	Prox. Humerus pl. EIGHT ² vas SL T IMPLANTS sterile
19.3683.10	Prox. Humerus plate EIGHT ² vas Corti SL T SET
19.3683.10PS+	Prox. Humerus pl. EIGHT ² vas Corti SL T IMPLANTS sterile
19.3751.10	VaWiKo [®] Humerus double threaded 3.5/3.5 Corti SL SET
19.3751.10PS+	VaWiKo [®] Humerus double threaded 3.5/3.5 Corti SL IMPLANTS sterile
19.3752.10	VaWiKo [®] Humerus double threaded 4.0/3.5 Corti SL SET
19.3752.10PS+	VaWiKo [®] Humerus double threaded 4.0/3.5 Corti SL IMPLANTS sterile
19.3753.10	VaWiKo [®] Humerus double threaded 4.0/3.5 Cancell SL SET
19.3753.10PS+	VaWiKo [®] Humerus double threaded 4.0/3.5 Cancell SL IMPLANTS sterile



★ Six lobe
Set no. 19.035.30I
Level 1
non-sterile packaging



★ Six lobe
Set no. 19.035.30I
Level 2
non-sterile packaging